
Mhat Ho 105 355 311

Week 2: Oscilloscopes and Function Generators

RMS of a periodic signal is calculated by first squaring the waveform, then taking its mean over its period, T, then taking the square root. Its definition using the calculus is

RMS =
$$\sqrt{\frac{1}{T} \int_0^T f^2(t) dt}$$

As an example, we will derive the equation for RMS/Vpp for a sine wave. You will be asked to derive the equation for square waves and triangular waves in the pre-lab.



First,
$$f(t) = a \sin\left(\frac{2\pi t}{T}\right) \to RMS = \sqrt{\frac{1}{T} \int_0^T a^2 \sin^2\left(\frac{2\pi t}{T}\right) dt}$$

Using the definition of
$$\sin^2(\theta) = \frac{1-\cos(2\theta)}{2}$$
, $RMS = \sqrt{\frac{a^2}{T} \int_0^T \frac{1-\cos(\frac{4\pi t}{T})}{2} dt}$

Taking the integral,
$$RMS = \sqrt{\frac{a^2}{T} \left[\frac{1}{2} t - \frac{\frac{1}{2}T}{4\pi} \sin\left(\frac{4\pi t}{T}\right) \right]_0^T}$$

Evaluating, we get
$$RMS = \sqrt{\frac{a^2}{T} \left[\frac{1}{2}T\right]}$$
 (Note that at $t = 0$ and T, $\sin\left(\frac{4\pi t}{T}\right) = 0$)

Therefore,
$$RMS = \frac{a}{\sqrt{2}}$$
, and since $Vpp = 2a$, then $\frac{RMS}{Vpp} = \frac{1}{2\sqrt{2}}$


It may for the purposes of your lab helpful to think of RMS in terms of Vpp, like so:

$$RMS = \frac{Vpp}{2\sqrt{2}}$$

Week 2 Prelab


Calculate the ratio RMS/Vpp for the following signals. Show all your work!


Name:




1. Square Wave: RMS / Vpp = ?

a for
$$0 \le + < \frac{7}{2}$$








$$\int_{T}^{T} \int_{0}^{2} a^{2} dt = \int_{T}^{T} \int_{0}^{2} a^{2} \left(T-\delta\right) = 0$$

UID: 105355311

2. Triangular Wave: RMS / Vpp = ? 2(3)




3. If you see a difference by a factor of 10 between the oscilloscope reading and the function generator setting, where is the first place that you should look? Watch the Probe Setting video (https://youtu.be/dtSuTHIviSo) for the answer.

Hirstly we should check the probe attenuation by looking at the channel 1 (yellow) satio on the oscilloscope If the satio is 10:1, it means We are connecting a 10-to-1 attenuation probe to the oxcill sope that causes the wrong measurement. So, we need to fix it by 41 adjusting the know to have a channel 1's ration is 1:1

$$= | ba^{2} \int_{37/4}^{7} (\frac{t}{7} - 1)^{2} dt \quad \text{Let } x = \frac{t}{7} - 1$$

$$= | ba^{2} \int_{37/4}^{7} = | x^{2} - 1/4$$

$$= | ba^{2} \int_{-1/4}^{8} = | ba^{$$



5. Why would you ever want to use AC coupling on an oscilloscope? Watch the AC Coupling video (https://www.youtube.com/watch?v=dtSuTHIviSo&t=6s) for the answer.

Because the AC coupling can fifter out the DC Signal component from a AC-DC signal. This will increase the solution of the Signal measurement on the oscilloscope.

Week 2 Prelab End